

=> d his

(FILE 'HOME' ENTERED AT 12:38:41 ON 16 MAR 2009)

FILE 'CA' ENTERED AT 12:38:51 ON 16 MAR 2009

L1 24735 S PROTEIN(4A) (REMOV? OR FILTRAT? OR ULTRAFILTRAT? OR MICROFILTRAT?  
OR ULTRAMICROFILTRAT?) OR DEPROTEINAT?L2 10476 S L1 AND(FILTRAT?  
OR ULTRAFILTRAT? OR MICROFILTRAT? OR ULTRAMICROFILTRAT? OR FILTER?  
OR ULTRAFILTER? OR MICROFILTER? OR ULTRAMICROFILTER?)  
L3 54 S L2 AND (ACETONITRILE OR CH3CN)  
L4 30 S L3 AND PY<2003  
L5 11 S L3 NOT L4 AND PATENT/DT  
L6 41 S L4-5

=> d bib,ab 16 1-41

L6 ANSWER 14 OF 41 CA COPYRIGHT 2009 ACS on STN  
AN 137:43736 CA  
TI The comparison of plasma deproteinization methods for the detection of  
low-molecular-weight metabolites by 1H nuclear magnetic resonance  
spectroscopy  
AU Daykin, Clare A.; Foxall, Peta J. D.; Connor, Susan C.; Lindon, John C.;  
Nicholson, Jeremy K.  
CS Biological Chemistry, Biomedical Sciences Division, Faculty of Medicine,  
Imperial College of Science, Technology and Medicine, London, SW7 2AZ,  
UK  
SO Analytical Biochemistry (2002), 304(2), 220-230  
AB Blood plasma is the major vehicle by which metabolites are transported  
around the body in mammalian species, and chem. anal. of plasma can  
provide a wealth of information relating to the biochem. status of an  
individual and is important for diagnostic purposes. However, plasma is  
very complex in physicochem. terms because it is composed of a range of  
org. and inorg. constituents with a wide range of mol. wts. and chem.  
classes and this makes anal. non-trivial. It is now well established  
that high-resoln. 1H NMR spectroscopy of blood plasma provides useful  
qual. and quant. biochem. information relating to metabolic disorders.  
However, one of the problems encountered in NMR spectroscopic anal. of  
blood plasma is the extensive peak overlap or presence of broad  
macromol. peaks in the 1H NMR spectrum, which can severely limit the  
amt. of obtainable information. Even with spectroscopic editing,  
information relating to low-mol.-wt. (MW) metabolites is frequently  
lost. Therefore, the efficiency of a range of conventional **protein  
removal** methods, in combination with the use of one- and two-dimensional  
NMR spectroscopic methods for evaluation, have been compared for the  
extn. of NMR-observable low-MW metabolites. It has been shown that  
these "deproteinization" methods vary considerably in recovery of low MW  
metabolites and a judicious choice is crucial for optimal extn. of a  
given analyte. The results presented here show that while  
**ultrafiltration** provides the "safest" method of plasma deproteinization,  
the signal-to-noise ratio of the resultant 1H NMR spectra is poor. On  
the other hand, **acetonitrile** pptn. at physiol. pH allows the detection  
of more low-MW metabolites and at higher concns. than any other method  
and provides the further advantages of being a rapid and simple  
procedure.

L6 ANSWER 17 OF 41 CA COPYRIGHT 2009 ACS on STN  
AN 132:87674 CA  
TI Automated **protein** precipitation by **filtration** in the 96-well format  
AU Biddlecombe, R. A.; Pleasance, S.  
CS Department of International Bioanalysis, Division of Bioanalysis and  
Drug Metabolism, GlaxoWellcome R&D, Ware, UK  
SO Journal of Chromatography, B: Biomedical Sciences and Applications  
(1999), 734(2), 257-265  
AB The use of automated **protein** pptn. by **filtration** in the 96-well format  
as a rapid sample prepn. technique for high throughput bioanal. using  
liq. chromatog. tandem mass spectrometry is reported. A robotic sample  
processor is used to aspirate sequentially a plasma sample and  
**acetonitrile** sepd. by air gaps. These are then mixed by being dispensed  
into individual channels of a 96-well **filter** block. The resulting  
supernatant is sepd. from the pptd. plasma proteins by the application  
of gentle vacuum using a custom manifold. The **filtered** supernatants are  
collected into a deep well microtiter plate, evapd. to dryness using a  
heated 96-well dry down station and reconstituted in water prior to  
anal. The efficiency of the extn. procedure is measured by the Lowry  
method for detg. protein concn. This method was used to optimize both  
the vol. and the order of reagent addn., and to compare several  
prototype 96-well **filter** blocks. Using the optimized procedure a  
specific, precise and accurate method was developed for the  $\beta$ -agonist  
salbutamol in rabbit plasma with a calibration range of 1 to 100 ng/mL  
from 100  $\mu$ l of sample.

L6 ANSWER 20 OF 41 CA COPYRIGHT 2009 ACS on STN  
AN 130:90010 CA  
TI Therapeutic drug monitoring of doxorubicin in pediatric oncology using  
capillary electrophoresis  
AU Hempel, Georg; Schulze-Westhoff, Petra; Flege, Silke; Laubrock, Nora;  
Boos, Joachim  
CS Kinderklinik, Abteilung Haematologie/Onkologie, Univ. Muenster,  
Muenster, D-48129, Germany  
SO Electrophoresis (1998), 19(16-17), 2939-2943  
AB A method for the detn. of doxorubicin and its main metabolite  
doxorubicinol in blood plasma is described. 2 Different sample prepn.  
procedures are applied depending on the expected concn. To monitor the  
peak plasma levels, 10  $\mu$ L of plasma are **deproteinated** with **acetonitrile**.  
After centrifugation, the supernatant is directly applied to the  
capillary by hydrodynamic injection. For the detn. of lower amts. of  
doxorubicin and doxorubicinol 100  $\mu$ L of plasma is extd. by liq.-/liq.  
extn. with chloroform. After evapn. of the org. phase, the sample is  
reconstituted in **acetonitrile**/water (95/5 vol./vol.) and injected into  
the capillary by electrokinetic injection. Idarubicin serves as the  
internal std. Laser-induced fluorescence detection with an Ar-ion laser  
emitting at 488 nm and a 520 nm cut-off **filter** is used for detection.  
The accuracy of the method was calcd. to be 3.0% at higher concns. and  
15.0% at the limit of quantification. Reproducibility data are in  
accordance to the generally accepted criteria for bioanal. methods. The

limit of quantification is 2 µg/L, enabling us to monitor doxorubicin plasma levels for several days after application. Noninvasive blood sampling (from the fingertip) using heparinized capillaries was a simple and convenient procedure and provides reproducible data. Initial results show high interindividual variability in doxorubicin peak plasma levels.

L6 ANSWER 22 OF 41 CA COPYRIGHT 2009 ACS on STN

AN 126:248477 CA

OREF 126:47983a,47986a

TI A study of deproteinization methods for subsequent serum analysis with capillary electrophoresis

AU Ralston, Pamela B.; Strein, Timothy G.

CS Dep. Chem., Bucknell Univ., Lewisburg, PA, 17837, USA

SO Microchemical Journal (1997), 55(2), 270-283

AB Methods of serum **protein removal** are examd. for applicability with subsequent serum anal. for low mol. mass metabolites and pharmaceutical agents by capillary electrophoresis (CE). Each of the deproteinization methods considered here has been reported in the literature and successfully used for serum **protein removal** in other applications. Two of the methods are found to be superior to the others for CE anal. **Protein removal** by **ultrafiltration** with Centrifree Micropartition System **filters** and by pptn. using **acetonitrile** worked well for **protein removal**, as evidenced by reproducible (RSD < 6%) detns. of several low mol. wt. components. The former is capable of analyzing for free metabolites only, while the latter was found to be somewhat problematic with respect to evapn. of the **acetonitrile**, causing potential problems with reproducibility. Both techniques for **protein removal** were clearly suitable for the subsequent serum anal. by CE in untreated fused silica capillaries.

L6 ANSWER 24 OF 41 CA COPYRIGHT 2009 ACS on STN

AN 125:8898 CA

OREF 125:2023a,2026a

TI Determination of pantothenic acid in infant milk formulas by high performance liquid chromatography

AU Romera, J. M.; Ramirez, M.; Gil, A.

CS Res. Dev. Dep., Abbott Laboratories S.A., Granada, 18004, Spain

SO Journal of Dairy Science (1996), 79(4), 523-526

AB A reverse-phase liq. chromatog. method was adapted for the assay of pantothenic acid in infant milk formulas. Sample prepn. consisted of **deproteination** with acetic acid and sodium acetate solns., followed by centrifugation and **filtration**. The chromatog. system included a C-18 column and a mobile phase consisting of a sodium phosphate buffer and **acetonitrile** (97:3, vol/vol). The column effluent was monitored by UV detection at 197 nm. The system was linear from 50 to 800 ng. The recoveries of pantothenic acid from augmented samples ranged from 89 to 98%, and the coeffs. of variation ranged from 1.17 to 3.20%. The results obtained with the HPLC and a microbiol. method were highly correlated for starting infant formula, followup infant formula, and formula for infants of low birth wt. from four different manufacturers. All formulas analyzed contained pantothenic acid at concns. higher than

those declared on their nutritional labels and were in compliance with international recommendations.

L6 ANSWER 32 OF 41 CA COPYRIGHT 2009 ACS on STN  
AN 113:90812 CA  
OREF 113:15087a,15090a  
TI Measurement of verapamil in human plasma by reversed-phase high-performance liquid chromatography using a short octyl column  
AU Rustum, Abu M.  
CS Dep. Environ. Fate Metab., Hazleton Lab. America, Inc., Madison, WI, 53707, USA  
SO Journal of Chromatography, Biomedical Applications (1990), 528(2), 480-6  
AB An HPLC method for the detn. of the Ca blocker verapamil in human blood plasma is described. The samples were **deproteinated** with **acetonitrile**/ZnSO<sub>4</sub>/MgSO<sub>4</sub>, centrifuged, and **filtered** prior to the anal. on a reversed-phase C-8 column. The anal. used a mobile phase of **acetonitrile**-phosphate buffer (60:40) adjusted to pH 7.1. The UV-detector was set at 220 nm. The calibration was linear in the range of 20-500 ng/mL, the accuracy was 95.0-98.3%, and the coeffs. of variation were 2.4-5.1%. No interference from plasma components and common drugs was found.

L6 ANSWER 33 OF 41 CA COPYRIGHT 2009 ACS on STN  
AN 111:70204 CA  
OREF 111:11655a,11658a  
TI A high-performance liquid chromatographic method for the measurement of deferoxamine in body fluids  
AU Tesoro, Angelo; Leeder, J. Steven; Bentur, Yedidia; Klein, Julia; Freedman, Melvin; Koren, Gideon  
CS Res. Inst., Hosp. Sick Child., Toronto, ON, M5G 1X8, Can.  
SO Therapeutic Drug Monitoring (1989), 11(4), 463-70  
AB A high-performance liq. chromatog. method for the anal. of deferoxamine (DFO) in 100 µL of serum or plasma is described. The procedure involves the addn. of the internal std. ciprofloxacin to the sample, followed by **ultrafiltration** to **remove protein**. The **ultrafiltrate** is then directly injected into the chromatog. system. Sepn. is achieved using a reverse-phase µBondapak C18 column and a ternary solvent system (sodium phosphate:**acetonitrile**:methanol) running at 2.0 mL/min. Assay time is 10 min, and chromatograms show no interference from coadministered drugs during this period of time. Coeffs. of variation were found to be less than 5%, and anal. recovery of DFO was 85%. Validation expts. in an exptl. dog model and in patients with iron overload demonstrate that the method is appropriate for studying the pharmacokinetics of DFO in thalassemic patients receiving drug for the treatment of chronic iron overload.

L6 ANSWER 38 OF 41 CA COPYRIGHT 2009 ACS on STN  
AN 107:112183 CA  
OREF 107:18127a,18130a  
TI Deproteinizing methods evaluated for determination of uric acid in serum by reversed-phase liquid chromatography with ultraviolet detection  
AU Sakuma, Ryoza; Nishina, Toshihiro; Kitamura, Motoshi

CS Dep. Clin. Chem., Toranomon Hosp., Tokyo, 105, Japan  
SO Clinical Chemistry (Washington, DC, United States) (**1987**), 33(8), 1427-30  
AB Six deproteinizing methods were evaluated for detn. of uric acid in serum by high-performance liq. chromatog. with UV detection: those involving zinc hydroxide, sodium tungstate, trichloroacetic acid, perchloric acid, **acetonitrile** and centrifugal **ultrafiltration** (with Amicon MPS-1 devices). A Toyosoda ODS-120A reversed-phase column was used. The mobile phase was sodium phosphate buffer (40 mmol/L, pH 2.2) contg. 20 mL of methanol per L. Absorbance of the eluate was monitored at 284 nm. The pptn. method with perchloric acid gave high recoveries of uric acid and good precision, and results agreed with those by the uricase-catalase method of N. Kageyama (1971).

=> log y

STN INTERNATIONAL LOGOFF AT 12:47:20 ON 16 MAR 2009